

Endocytosis of Siglec-8 and Siglec-F on bone marrow-derived eosinophils from a ROSA26-Siglec-8-KI EPX-Cre mouse

- Jeremy O'Sullivan

Materials:

- Bone marrow-derived eosinophils (BMDEs) at day 10-14 of *in vitro* differentiation
- Flow buffer (autoMACS Running Buffer; Miltenyi Biotec, cat# 130-091-221)
- Unconjugated anti-Siglec-8 (2C4)
- Unconjugated anti-Siglec-F (in house hybridomas)
- FITC-labeled anti-mouse IgG1 secondary Ab (BioLegend, cat# 406606)
- PE-labeled anti-rat IgG secondary Ab (BD Pharmingen, cat# 550767)
- Fluoropure-grade DAPI (as viability stain; ThermoFisher Scientific, cat# D21490)
- 5-ml polystyrene round-bottom flow tubes (Falcon, ref. no. 352052)
- RPMI-1640 medium (Gibco, cat# 72400-047)
- BD LSR II flow cytometer
- FlowJo analysis software


Experimental Procedure:

1. Transfer 2×10^5 to 1×10^6 total cells to each of the flow tubes for staining. Separate tubes will be needed for:
 - a. An unstained sample
 - b. A sample stained only for Siglec-8 (single stain control for compensation)
 - c. A sample stained only for Siglec-F (single stain control for compensation)
 - d. A sample stained only with DAPI (single stain control for compensation)
 - e. An isotype control sample with isotype controls for both anti-siglec antibodies to identify non-specific and background staining
 - f. One sample for each time point of endocytosis, including one at 0 minutes (to indicate initial surface levels of Siglec-8 and Siglec-F)
2. Wash with 1-2 ml of flow buffer and spin down the cells at $300 \times g$ (1400 RPM) for 5 min at 4°C.
3. Decant the supernatant and vortex to resuspend the cells.

4. **Unconjugated primary antibody labeling:** Incubate the samples in 100 μ l of flow buffer with the appropriate antibody/antibodies for 20 min at 4°C. To do this, make stock 2X dilutions of each of the unconjugated antibodies (5 μ g/ml of the anti-Siglec-8 mAb [2C4] and 10 μ g/ml of the anti-Siglec-F mAb, the same concentrations for the isotype controls) in flow buffer with 50 μ l for each applicable sample, plus about 10% for pipetting error. Add 50 μ l of each dilution to the relevant sample tubes and bring to 100 μ l with flow buffer, if necessary.
5. Following the 20-min incubation, wash with 1-2 ml of flow buffer and spin down the cells at 300 x g (1400 RPM) for 5 min at 4°C.
6. Decant the supernatant, resuspend the cells, and add 100 μ l of cold RPMI-1640.
7. **Endocytosis incubation step:** Keep the controls and 0-min sample on ice and incubate the experimental samples for the appropriate duration in the incubator at 37°C.
8. Wash all of the samples with 1-2 ml of flow buffer and spin down the cells at 300 x g (1400 RPM) for 5 min at 4°C. While the samples are spinning, make up dilutions of the detection antibodies. Make 2X dilutions of the FITC-conjugated anti-mouse IgG1 (1:100, for a final dilution of 1:200; to detect the anti-Siglec-8 or its isotype control) and the PE-conjugated anti-rat IgG (1:100, for a final dilution of 1:200; to detect the anti-Siglec-F or its isotype control) with 50 μ l for each applicable sample, plus about 10% for pipetting error.
9. **Fluorophore-conjugated secondary antibody labeling:** Decant the supernatant, resuspend the cells, and add 50 μ l of each dilution to the relevant sample tubes and bring to 100 μ l with flow buffer, if necessary. Incubate for 20 min at 4°C.
10. Wash all of the samples with 1-2 ml of flow buffer and spin down the cells at 300 x g (1400 RPM) for 5 min at 4°C.
11. Decant the supernatant, resuspend the cells, and add 200 μ l of DAPI-containing flow buffer (pre-diluted, at 1.8 μ M).
12. Analyze the samples by flow cytometry.


Endocytosis of Siglec-8

...on mouse Siglec-8+ bone marrow-derived eosinophils over 4 hours:




Endocytosis of Siglec-8

...on mouse Siglec-8+ bone marrow-derived eosinophils over 4 hours:




Endocytosis of Siglec-8

...on mouse Siglec-8+ bone marrow-derived eosinophils over 4 hours:




Endocytosis of Siglec-8

...on mouse Siglec-8+ bone marrow-derived eosinophils over 4 hours:




Endocytosis of Siglec-8

...on mouse Siglec-8+ bone marrow-derived eosinophils over 4 hours:




Endocytosis of Siglec-8

...on mouse Siglec-8+ bone marrow-derived eosinophils over 4 hours:




Endocytosis of Siglec-8

...on mouse Siglec-8+ bone marrow-derived eosinophils over 4 hours:



Endocytosis of Siglec-8

...on mouse Siglec-8+ bone marrow-derived eosinophils over 4 hours:



240 min