

Module 3: Preparation of liposomes and their use in flow cytometry

Objective:

Make human CD22 targeted and Naked liposomes and to test their binding to human CD22 expressing CHO cells using flow cytometry.

References:

Chen WC et al, In vivo targeting of B-cell lymphoma with glycan ligands of CD22. *Blood*. 2010, 115(23): 4778–4786.

Reagents and Equipments:

- Lipid components
 - 1,2-distearoyl-*sn*-glycero-3-phosphocholine (DSPC, Avanti #850365)
 - Cholesterol (Sigma-Aldrich #C8667)
 - N-(Carbonyl-methoxypolyethyleneglycol 2000)-1,2-distearoyl-*sn*-glycero-3-phosphoethanolamine (PEG-DSPE, NOF # SUNBRIGHT DSPE-020CN)
 - 1,2-dioleoyl-*sn*-glycero-3-phosphoethanolamine-N-(7-nitro-2-1,3-benzoxadiazol-4-yl) (ammonium salt) (NBD-PE, Avanti #810145C)
 - Human CD22 glycan ligand-PEG-DSPE
- Liposome extruder set (Avanti, #610000)
- Filter support (Avanti # 610014)
- Polycarbonate membranes (Avanti #610005 for 100 nm, 610007 for 400 nm)
- Human CD22-CHO and parental CHO cells
- HBSS + 0.01%BSA (FACS buffer)
- PBS without Ca and Mg
- Glass tube (Fisher #14-958D)
- FACS tube (BD Biosciences #352052)

Instruments:

- Sonication bath (BRANSON ultrasonic bath 1510-MTH)
- Flow cytometer (FACS caliber, BD Biosciences)
- (Optional) Particle size analyzer (Zeta sizer, Malvern)

Methods:

1. Liposome preparation:

- Preparation of lipid components (1mmol total lipids)

DSPC (MW=790)	Chol (MW=387)	PEG-DSPE (MW=2900)	CD22 ligand-PEG-DSPE (MW=3880)	NBD-PC (MW=856)
56mol%	38mol%	5-Xmol%	Xmol%	1mol%
0.56 mmol	0.38 mmol	0.05-0.0X mmol	0.0X mmol	0.01 mmol

- Dissolve the DSPC, Chol, PEG-DSPE, and NBD-PC in Chloroform (CD22 ligand-PEG-DSPE is in DMSO).
- Add lipid components (except CD22 ligand) in glass tube.
- Vortex the tube.
- Dry the chloroform with N₂ gas flow.
- Add CD22 ligand-PEG-DSPE in DMSO.
- Lyophilize the DMSO for 18 h.

Norihito will prepare the dried lipid components before the day.

- Hydration of lipid components
 - Add 1 mL of PBS into glass tubes (liposome concentration is 1 mM).
 - Sonicate the solution for 3 min with holding the glass tube.
 - Vortex the tube occasionally.
- Extrusion of liposomes
 - Set up the extruder according the following link.
http://avantilipids.com/index.php?option=com_content&view=article&id=531&Itemid=295
 - Pass the membrane filter gently over 20 times using 400 nm and 100 nm pore size filters.
 - (Optional) Measure the size of liposomes with the light scattering analyzer zetasizer.

2. Analysis of liposome binding to human CD22-CHO cells:

- Staining of the cells with liposomes
Cells (1-3 x 10⁵ cells in 100 mL) in FACS tube (**Norihito will prepare in advance**)

- Add 100 mL of liposome solution (**Dilution of liposome will be determined by Norihito in advance**).
 - Vortex the FACS tube gently.
 - Incubate the cells for 30 min at 37°C.
 - Wash the cells by adding 2 mL of FACS buffer (No need for pipetting).
 - Spin down the tubes at 510 x g for 3 min.
 - Discard the supernatants by pouring the supernatants into waste.
 - Vortex the tube to loosen the cell pellet.
 - Suspend the cells in 300 mL of FACS buffer containing 1 mg/mL propidium iodide and transfer the cells into the FACS tubes.
 - Vortex the tube well.
- Analyze the liposome binding by FACS Caliber.
 - Turn on the FACS caliber.
 - Turn on the computer.
 - Open the Cell quest software.
 - Make a dot plot of FSC vs SSC (This is for cell detection).
 - Make a dot plot of FSC vs FL3 (This is for dead cell exclusion).
 - Make a histogram plot of FL1 (This is for NBD detection).
 - Make a folder to save the acquired files.
 - Acquire each sample.

For further information of FACS Caliber, please download the tutorials from the link below.

<http://www.bdbiosciences.com/support/resources/facscalibur/index.jsp>

- Analyze the data with Flowjo software.
 - Open the Flowjo software.
 - Drag and drop the folder containing the acquired samples on the flowjo.
 - Open the file of “unstained cells” and make a gate to identify the cells on the FSC vs SSC plot.
 - Open the gated population and make a dot plot of FSC vs FL3.

- Make a new gate for the living cells (FL3 low).
- Open the gated population and make the histogram plot of FL-1.
- Drag the newly generated files on the “all samples”. Flowjo will automatically make the same gate and histogram plot for all the samples.
- Open the work space in flowjo to draw the histogram data.
- Drag and drop the FL-1 histogram of “unstained cells”.
- To make a overlay histogram, drag and drop the “naked” and “CD22-tageted” files on the histogram of “unstained cells”.

For further information of Flowjo, please download the tutorials from the linke below.

<http://www.flowjo.com/home/tutorials/>

- Make slides in PowerPoint.
 - Copy and paste the overlay histogram on powerpoint file.
 - Make a figure legend.